

Bakezyme-GO SUP



Bakery Enzymes

Glucose Oxidase

PRODUCT DESCRIPTION

Bakezyme-GO SUP is a glucose oxidase enzyme that strengthens gluten structure, enhancing dough handling and machinability.

It helps in providing better dough consistency, thereby making it easier to process and shape during production.

With the enhancement of the gluten network, this enzyme also contributes to better texture and rise in baked goods, thus leading to higher quality breads and other bakery products. This is an essential ingredient for bakers looking to achieve optimal dough performance and uniform product quality.

KEY FUNCTIONS

- Strengthens gluten structure.
- Improves dough handling and machinability.

TYPICAL APPLICATIONS

Pan bread & buns

Flatbreads & rotis

Cakes & sweet goods

Biscuits & crackers

Flour milling & correction

Frozen dough

HANDLING & DOSAGE

Optimum dosage depends on substrate, process conditions (pH, temperature, time) and the desired result. We recommend lab or pilot trials; our application team can suggest a starting dose. Handle as an enzyme/biological preparation: avoid dust formation and inhalation, and wear gloves and eye protection.

SPECIFICATION

Product	Bakezyme-GO SUP (Glucose Oxidase)
Product range	Bakery Enzymes
Function	Glucose Oxidase
Enzyme activity	Specified per grade — on request
Physical form	Powder; other formats on request
Packaging	25 / 50 kg HDPE drums or retail pouches
Storage	Cool, dry place in tightly closed original containers, away from direct sunlight
Documents	Certificate of Analysis (CoA) on request

MANUFACTURING & QUALITY

Manufactured by Enzyme Bioscience in an FDA-registered, GMP and ISO 9001:2015 certified facility in Surat, India, with in-house R&D and QC laboratories.



GMP

ISO
9001:2015FDA
registered

The information in this data sheet describes typical properties and is given in good faith; it does not constitute a specification or warranty for a particular use. Users should confirm suitability for their own process and comply with applicable regulations. Batch-specific values are stated on the Certificate of Analysis.